



*Be Right™*

# How and Why to Monitor ATP for Biological Growth

Ted Simmons – Regional Sales Manager, Southern Ohio  
tsimmons@hach.com

Company Confidential

# Agenda

- Overview of the Hach Company
- Hach and LuminUltra
- What is ATP
- Options for testing
- What is the value of knowing the microbial load
- Advance testing from LuminUltra
- Process ATP testing is EZ

# Overview of Hach Company

- **Began in 1947 in Ames, IA**
- **1999: Acquired by Danaher Corporation**
- **Manufacturer of industrial and municipal water analysis solutions**
  - process and laboratory instruments
  - chemistries
  - service and software
- **Innovation leader**
  - 527 patents, 130 patent families
  - strong investment in R&D and acquired technologies



- **OUR MISSION**

Ensure water quality for people around the world.


- **OUR VISION**

We make water analysis better—faster, simpler, greener and more informative—via unsurpassed customer partnerships, the most knowledgeable experts, and reliable, easy-to-use products.



# Overview of Hach Company

- **70+ RSM's**
  - Regional Sales Managers
- **25+ CAM - KAM – ADM – TSS - CSS**
  - CAM - Corporate Account Managers
  - KAM - Key Account Managers
  - ADM - Application Development Managers
  - TSS - Technical Sales Specialist
  - CSS - Complex Sales Specialist
- **100+ Field Service Technicians**
  - Field Service Partnerships
  - Startups, Commissioning and Training
- **41+ Technical Support Representatives**
  - Global Customer Support
- **12+ Bench Service Technicians**
  - Two service centers for repair and certification





# Hach and LuminUltra Teaming Up



## May of 2013

- Hach is a non-exclusive distributor of the LuminUltra Product Line
- Tech Support will be prepared to answer questions related to the following
  - Pre-Sale questions related to LuminUltra product offering
  - Ability to identify customers that may be interested in ATP testing
  - Ability to specify the correct LuminUltra test kits for different applications
  - Post-Sale troubleshooting of LuminUltra product offering
  - Exhibit elite level of support, ownership, and product knowledge LuminUltra product line



# About LuminUltra

**Founded in 1995**, LuminUltra is a biological diagnostic testing company headquartered in Canada with operations in 6 countries. It is widely recognized globally as a leader in developing tests and reagents for environmental, industrial, and diagnostic monitoring and is a key supplier of COVID-19 clinical testing reagents to the Government of Canada. Customers in over 80 countries trust LuminUltra's technology, production reliability and history of customer service excellence to deliver their essential services in a safe state.

LuminUltra proudly serves some of the top companies in the world including:



Baltimore

Fredericton

London

Paris

Heerenveen

Melbourne

6

LuminUltra offices spread across the globe

3.5K

Customers from around the world leverage LuminUltra's testing solutions

80+

Countries with customers that rely on LuminUltra's products and services

\$345M

In customer value delivered to date

25

Years of exceptional production reliability and innovation

500K


Clinical COVID-19 testing reagents per week being supplied to Canada

# The Water Quality Toolbox

- Many physical and chemical parameters can be determined in minutes (or instantly via online analyzers).
- The technology has not existed to allow microbiological analyses to be done in the field and produce results quickly.

© 2020 LuminUltra Technologies Ltd.

- This results in significant limitations in water quality control capabilities.



# The Problem

- Several Hach products give water testers the ability to quickly and easily assess several water quality parameters, such as:
  - Temperature
  - pH
  - Alkalinity
  - Turbidity
  - Color
  - TDS
- But what about quick assessments about biological activity? Options are **very limited.**






# What is ATP?

- ATP = Adenosine Triphosphate
  - Primary energy carrier for all life.
  - Measurement takes only minutes.
- Since ATP is present in ALL living cells, quantifying it enables you to detect the entire population rather than just culturable cells.
- ATP testing has been used in food & medical hygiene for decades.

Company Confidential







Be Right™

# The Basics of ATP Monitoring


- **ATP = Adenosine Triphosphate**
  - Primary energy carrier for all life.
  - Measurement takes only minutes.
- LuminUltra is the **first** and **only** company to provide 2<sup>nd</sup> Gen. technology designed for fluid systems (far more sensitive, robust and reliable).





# Industry Standard – HPC

- HPC's incubate for 24-48 hours.
- **HPC's detect only ~0.1-1% of bioburden.**
- Information is only provided on organisms that can grow...
  - ...in the media used;
  - ...at the temperature provided;
  - ...within the incubation time allowed.
- Regulatory parameters (e.g. Total Coliforms, E. Coli) measure a even smaller fraction.
- **Many problem microorganisms are missed!**




# Strategy: Traditional Culture Tests



© 2020 LuminUltra Technologies Ltd.

## Strategy: In-field ATP testing



Point-and-shoot strategy can be applied to isolate the root cause and solve problems on-the-spot.

# What are the results of uncontrolled microbial growth?

- Biofilm accumulation
- Elevated disinfectant demand
- Taste & odor, cosmetic issues
- Non-compliance
- Premature infrastructure degradation (i.e. MIC)

© 2020 LuminUltra Technologies Ltd.

***Which of these can be managed using compliance tests?***

Company Confidential

# Limitations of Traditional Microbiological Testing

- None of the problems discussed on the previous slides are traced to a single type of organism – they are caused by the community as a whole.
- Attempts at total microbial content (e.g. HPC) are slow and a poor representation of the total population.
- The tools being used are not appropriate for the objective!

© 2020 LuminUltra Technologies Ltd.

***“90-99.9% of the bacteria found in water systems are not culturable in common laboratory media.”***

Risenfeld, 2004





# What are the **costs** of uncontrolled microbial growth?

- Clean ≠ Compliant! Significant costs are incurred as a result of general growth.
- Several hours and large volumes of water required for flushing programs. Adds up to **>\$122 US/hour** to flush.

Assumptions:

- Water value = \$2 US/1000 USgal
- Hydrant flow rate = 75 USgpm
- **\$1/person/day** for boil water advisory (conservative estimate).

© 2020 LuminUltra Technologies Ltd.

Company Confidential

# Costs of Line Breaks

Location	Miles of Pipe	Total # of Breaks	Cost to Manage Annual Breaks <sup>1</sup>
Phoenix	5,400	1350	\$4,050,000
Dallas	4,600	1150	\$3,450,000
Kansas City	1,912	478	\$1,434,000
Macon	1,400	350	\$1,050,000

- **Average line break frequency:**  
23-27 breaks/100 miles of pipe/year (Kirmeyer et al)

© 2020 LuminUltra Technologies Ltd.

- **Direct costs to repair a break: \$3000**  
(American Water Works Service Company)

Company Confidential



# Costs of Line Breaks

Location	Miles of Pipe	Total # of Breaks	Cost to Manage Annual Breaks <sup>1</sup>	Potential Savings
Phoenix	5,400	1350	\$4,050,000	\$405,000
Dallas	4,600	1150	\$3,450,000	\$345,000
Kansas City	1,912	478	\$1,434,000	\$143,400
Macon	1,400	350	\$1,050,000	\$105,000

- Assume 50% of breaks due to corrosion (the balance is due to construction incidents, etc.)

© 2020 Lumina Ultra Tech, Inc.

- Also (conservatively) assume that enhanced control reduces field service time by 10%.

Company Confidential



## As a result...

- People have learned to live without microbiological data.
  - “We don’t do micro. That’s for the people across town/state lab/etc.”
  - “We just keep our chlorine high.”
  - “We just look at Total Coliform/E.Coli. We don’t really care about other stuff.”
- Many are aware of limitations but aren’t aware that a solution exists.

## As a result...

- If a problem is detected using traditional microbiological tools, it has been there for at least a couple of days.
- Since the time of sample collection, the problem may (and likely is) worse and spread downstream.
- What would it be worth to you to know about the problem immediately?



# The Solution



- **Value Proposition** – A rapid, non-specific measure of living microorganisms in water (or any other fluid sample).
- Four main advantages of LuminUltra's tests:
  1. Real-time feedback (< 5 minutes)
  2. Complete results (100% of species detected).
  3. Field-ready
  4. Test Kits designed for specific applications
- Decisions can be made on-the-spot, enabling **same-shift troubleshooting**.



## PRODUCT OVERVIEW


# 2<sup>nd</sup> Generation ATP

### “The first line of defence”


- Trusted solution for water-related industries for 15+ years;
- Rapid measure of total bioburden in < 5 minutes using a simple and portable test protocol;
- Facilitates timely risk assessments to identify areas of concern;
- Drives immediate action (i.e. more specific assays and/or adjustments to treatment strategies).



Non-Filterable Samples: cATP (Cellular ATP) = tATP (Total ATP) - dATP (Dissolved ATP)




Filterable Samples: cATP (Cellular ATP) Direct






- Luminometer-PhotonMaster
  - Powered via USB to laptop or desktop computer
- Application Specific Test Kits
  - Quench-Gone Aqueous (QGA)
  - Quench-Gone Organic Modified (QGO-M)
  - QuenchGone21 Wastewater (QG21W)
  - QuenchGone21 Industrial (QG21I)
  - QuenchGone21 Specialty (QG21S)
  - DSA Test Kit



- LumiCalc Software
  - Significant time savings in the organization, analysis, and reporting of data – spend minutes to achieve what took hours with traditional method.
  - Get interpretation as results are collected to make faster and better decision.



Company Confidential • Direct connection to PhotonMaster saves time tabulating results and protects your data

# Hach plus LuminUltra

**QuenchGone**  
The Next Generation of Microbial Monitoring  
**AQUEOUS**



**QuenchGone21**  
Microbial Monitoring for the 21<sup>st</sup> Century  
**WASTEWATER**



**QuenchGone21**  
Microbial Monitoring for the 21<sup>st</sup> Century  
**INDUSTRIAL**



**QuenchGone21**  
Microbial Monitoring for the 21<sup>st</sup> Century  
**SPECIALTY**



**Deposit&Surface**  
The Next Generation of Microbial Monitoring  
**ANALYSIS**



**QuenchGone**  
The Next Generation of Microbial Monitoring  
**ORGANIC-M**



Be Right™

Company Confidential

# 6 Test Kits + 1 Equipment Set

Be Right™

- ✓ Purified
- ✓ Drinking
- ✓ Cooling

**QuenchGone**  
AQUEOUS



**QuenchGone**  
ORGANIC-M



- ✓ Fuels
- ✓ Metalworking
- ✓ Oilfield

- ✓ Papermaking
- ✓ Pulping
- ✓ Wash Waters



- ✓ Polymers
- ✓ Slurries
- ✓ Adhesives



**QuenchGone21**  
SPECIALTY

- ✓ Deposits
- ✓ Biofilm
- ✓ Slimes
- ✓ Surfaces



- ✓ Wastewater
- ✓ Bioreactors
- ✓ Effluent
- ✓ Bioprocessing



### QGA Test Kit Contents & Storage Conditions

Component (LuminUltra P/N)	Storage	Shelf Life
<b>Luminase™ Enzyme &amp; Buffer Vials (Lu-3mL-FD)</b> <i>Luciferase Enzyme Reagent, 3mL</i>	4 to 25°C	6 to 12 mo*
<b>UltraCheck™ 1 Dropper Bottle (UC1-5mL)</b> <i>1 ng ATP/mL Standard, 5mL</i>	4 to 25°C	18 mo
<b>UltraLyse™ 7 Bottle (UL7-125mL)</b> <i>ATP Extraction Reagent, 125mL</i>	4 to 25°C	18 mo
<b>UltraLute™ (Dilution) Tube, 9mL (ULu-9mL-50R)</b> <i>ATP Dilution Reagent, 9mL</i>	4 to 25°C	18 mo
Quench-Gone Syringe Filters, 25/pk (DIS-SFQG-25)	-	-
60mL Syringe, PP/Neoprene, 25/pk (DIS-S60-25)	-	-
100 to 1000µL Blue Pipet Tips, 100/rack (DIS-PT1-100R)	-	-
10 to 200µL Yellow Pipet Tips, 96/rack (DIS-PT01-96R)	-	-
12x55mm Test Tubes, 50/pk (DIS-CT12-50)	-	-

**QGA** – Quench-Gone Aqueous. Use for water samples whether fresh or brackish from any source when microbial control is required.

Company Confidential

# QGA Quick-Reference Guide

## Step 2 - Cellular ATP (cATP™) Analysis

### 2.1 – MEASURE SAMPLE VOLUME

Determine volume and filter sample.

Sample Type	Volume (mL)
Cooling or Process Water	10 to 25
Fresh Brackish & Salt Water	25 to 50
Reclaimed Water, Effluents	25 to 50
Drinking and Sanitary Water	50 to 100
High Purity Water	> 100




### 2.2 – FILTRATION

Filter sample.



### 2.3 – EXTRACTION

Extract ATP from filter & dilute.



### 2.4 – ASSAY

Measure ATP concentration.



NOTE: If  $RLU_{cATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low-detection limit

NOTE: If  $RLU_{cATP} \leq 50$  using a PhotonMaster or Lumitester C-110, consider accounting for background ( $RLU_{bg}$ ). See Test Kit Instructions for guidance.



# QG21W Test Kit



**QG21W Test Kit Contents & Storage Conditions**

Component (Part Number)	Store At	Shelf Life
Luminase™ <sup>W</sup> Enzyme & Buffer Vials (LuW-3mL-FD) Luciferase Enzyme Reagent, 3mL	4 to 25°C	6 to 12 mo*
UltraCheck™ <sup>1</sup> Dropper Bottle (UC1-5mL) 1 ng ATP/mL Standard, 5mL	4 to 25°C	18 mo
UltraLyse™ <sup>3021</sup> (Extraction) Tube, 2mL (UL30(21)-2mL-50R) tATP Extraction Reagent, 2mL	4 to 25°C	18 mo
UltraLute™/Resin (Dilution) Tube, 8mL (ULuR-8mL-50R) tATP Dilution Reagent, 8mL	4 to 25°C	18 mo
LumiSolve™ (Stabilizer) Tube, 10mL (LS-10mL-50R) dATP Stabilizing Reagent, 10mL	4 to 25°C	18 mo
UltraLyse™ <sup>3021</sup> (Extraction) Tube, 10mL (UL30(21)-10mL-25R) ** tATP Extraction Reagent, 10mL	4 to 25°C	18 mo
100 to 1000µL Blue Pipet Tips, 100/rack (DIS-PT1-100R)	-	-
100 to 1000µL Wide-Mouth Pipet Tips, 100/rack (DIS-PT1WM-100R)	-	-
10 to 200µL Yellow Pipet Tips, 96/rack (DIS-PT01-96R)	-	-
12x55mm Assay Tubes, 50/pk (DIS-CT12-50)	-	-
2" x 2" 250µm Mesh Squares, 25/pk ** (DIS-MESH-25)	-	-
1.5" x 1.5" Weigh Boat, 25/pk ** (DIS-WD-25)	-	-
Scissor-Type Forceps, 1/pk** (EQP-FOR)	-	-

**QG21W** – QuenchGone21 Wastewater. Designed for quantification of biomass in wastewater treatment bioreactors, influents, and effluents. Also used for fermentation.

Company Confidential

**NOTE:** Kit uses **Luminase<sup>W</sup>**


**Luminase**



# QG21W Quick-Reference Guide

## Step 1 - UltraCheck™ 1 Calibration

Perform one UltraCheck 1 calibration per day or per each set of samples analyzed.




NOTE: If  $RLU_{ATP} \leq 500$  using a PhotonMaster or Lumitester C-110, rehydrate a new bottle of Luminase<sup>®</sup> for maximum sensitivity.

## Step 2 – Total ATP (tATP™)

Included in QG21W™ and QG21Wa™ test kits.


### 2.1 - EXTRACTION

Add sample to extract ATP.



### 2.2 – DILUTION

Dilute out interferences.



### 2.3 – ASSAY

Measure ATP concentration.




NOTE: If  $RLU_{tATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low- detection limit.

## Step 3 – Dissolved ATP (dATP™)

Included in QG21W™ and QG21Wa™ test kits.

### 3.1 – DILUTION

Add sample to recover ATP.



### 3.2 – ASSAY

Measure ATP concentration.




NOTE: If  $RLU_{dATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low- detection limit.

## STEP 4a – Floc Bulking ATP (fbATP™)

Included in QG21Wa™ test kit only.


### 4a.1 - FILTRATION

Filter sample to separate bulking floc.




### 4a.2 - EXTRACTION

Extract ATP from mesh.



### 4a.3 - DILUTION

Dilute out interferences.



### 4a.4 - ASSAY

Measure ATP concentrations.




NOTE: If  $RLU_{fbATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low-detection limit.

## STEP 4b – Attached Growth ATP (agATP™)


Included in QG21Wa™ test kit only.

### 4a.1 - EXTRACTION

Extract ATP from sample.




OR



### 4b.2 - DILUTION

Dilute out interferences.



### 4b.3 - ASSAY

Measure ATP concentrations.



NOTE: If  $RLU_{agATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low-detection limit.

*Municipal Water Management:*

## Actions driven by ATP test results

- Field-deployable test format enables user to trace up the line to find the source of issues;
- Optimize flushing time/water usage;
- Boost disinfectant/burnout;
- Conduct system audits to optimize secondary microbiological monitoring tools to speciate.





# Focus: Storage Tank Surveillance

- Storage tanks are often the first point at which regrowth becomes a problem.
- Whether it be due to long water age, stagnation, or infiltration, stored water represents a threat to downstream water quality.


$[ATP]_{outlet} > [ATP]_{inlet}$ ?

- If ATP leaving the tank is greater than that entering the tank then growth is occurring in the tank.




# Focus: Storage Tank Surveillance

- Water quality issues had been apparent through a water delivery line with the source suspected to be a storage tank on-site.
- Samples were drawn at the source, the tank inlet, and the tank outlet.
- ATP tests were done at each location over a 3-day period.



# Focus: Storage Tank Surveillance



© 2020 LuminUltra Technologies Ltd.

# Focus: Distribution System

1. Diagnose water quality issues as biological or not in minutes;
2. Trace the line to the source;
3. Apply corrective action;
4. Validate actions on-the-spot.

© 2020 LuminUltra Technologies Ltd.



Company Confidential

# Case Study: Distribution Troubleshooting

- A utility performed a distribution system audit to compare microbiological contamination levels with total chlorine residual.
- To recap: ATP targets for drinking water systems are as follows:

© 2020 LuminUltra Technologies Ltd.

**High Risk: >10pg/mL**


**Moderate:  $1 < \text{ATP} < 10$**

**Low Risk: < 1pg/mL**

Company Confidential

# Case Study: Distribution Troubleshooting

- As expected, ATP increased and total chlorine levels decreased with distance from the plant.
- BUT, in some cases ATP results persisted at elevated levels even in areas of adequately high chlorine.




© 2020 LuminUltra Technologies Ltd.

Company Confidential


# Case Study: Distribution Troubleshooting

- Chlorine was boosted in addition to a flushing cycle.
- Chlorine residuals increased as expected with the highest levels being nearest to the plant.
- After the flush, same levels of residual Cl were observed at all points.



# Case Study: Distribution Troubleshooting

- Despite the re-establishment of adequately high Chlorine, high ATP was still measured at the far ends of the distribution system.
- Localized focus could be placed in this area, whether it be in the form of additional flushing or additional root cause analysis.






# Case Study: Flushing

ATP levels dropped significantly after only 30 minutes to below the target 1 pg ATP/ml.

***Savings: 2 man hours +  
45,000 USgal of water.***

© 2020 LuminUltra Technologies Ltd.




# Case Study: Flushing

In this case ATP confirmed that the high turbidity was a non-biological issue and flagged the end of the flushing cycle after 30 minutes.

***Savings: 3 man hours +  
67,500 USgal of water.***

© 2020 LuminUltra Technologies Ltd.



# Focus: Nitrification

- Controlling nitrification requires a delicate balancing act:
  - $\text{NH}_3/\text{Cl}/\text{pH}$  – all a function of what's in the water, difficult to make sure it's all as it should be.
  - Nitrifying bacteria cannot be detected using HPCs.
  - HPCs detect Denitrifying bacteria, but by that time, it's already too late.



# Summary on ATP Testing

- ATP technology overcomes the limitations of traditional microbiology by taking microbiological testing into the field and producing results as quickly as other water quality parameters;
- Reveal, isolate and solve water quality problems in a fraction of the time that would otherwise be needed;
- Use results of in-field ATP tests to drive and prioritize additional testing to further characterize water quality.

## PRODUCT OVERVIEW

# GeneCount™ qPCR

### Isolates specific microbes

- qPCR testing has historically been complex, expensive, and a significant investment that only research labs would pursue.
- GeneCount is very “on brand” with LuminUltra’s traditional approach: practical and portable yet refusing to sacrifice accuracy for the sake of simplicity.
- Value Proposition: Test for the specific threats that concern you most without involving a lab, thereby saving significant time and money.
- Expandable platform to include numerous other tests in the future, including infectious diseases in addition to more pathogens and troublesome microbes.



## PRODUCT OVERVIEW

# GeneCount™ qPCR

### Additional Assays:

- Total Prokaryote
- Sulfate-Reducing Prokaryote
- Sulfur-Oxidizing Prokaryote
- Methanogens
- Iron Reducing Bacteria
- Total Fungi
- Total E.Coli
- Legionella sp.
- **NEW: SARS-CoV-2**


### Coming Soon:

- Legionella pneumophila
- Nitrifying Prokaryotes
- Nitrate-Reducing Bacteria
- Acid-Producing Bacteria
- Pseudomonas



Be Right™

# Having a microbial safety plan is crucial to ongoing risk management



LuminUltra's multifunctional GeneCount line provides both the platform and assays needed for ongoing pathogen defense!

# The EZ Series: online water analysis made easy

**ORGANIC LOAD ANALYZERS**

Total Organic Carbon  
Chemical Oxygen Demand

**NUTRIENT ANALYZERS**

Nitrate, phosphate, ammonia  
Total N, Total P

**TRACE METAL ANALYZERS**

12 heavy / trace metals  
Ag, As, Cr, Hg, Pb, Se

**NICHE WASTEWATER ANALYZERS**

Volatile fatty acids, FOS/TAC  
Influent Toxicity

**Matrix: WATER**



**MICROBIOLOGY ANALYZERS**

Adenosine Triphosphate (ATP)

**GENERAL WATER ANALYZERS**

Chemical parameters  
From A tot Z


**EZ SERIES:**  
Hach's advanced analyzer platform for online water analysis.

ANALYZER PORTFOLIO



# A Dynamic Process – Real Time Monitoring Needed

\* Graph by Montana State University



Planktonic / Sessile / Planktonic

## The Objectives:

- “Healthy biofilms will release lumps of bacteria once fully grown: sessile bacteria turn planktonic again: this is a risk.
- When Biocides are dosed they will only touch the surface of the biofilm releasing similar lumps of bacteria: risk.
- A clean system can be obtained by: cleaning and monitoring of low concentrations planktonic bacteria.


Cell dies and ATP is rapidly degraded

ATP serves as a reliable biomarker for living organisms in a sample


ATP can be used as a parameter in high-frequency, automatic, on-line microbial analysis in water

## Things to know about measuring ATP

Online options to measure Free ATP, Total ATP and Intracellular ATP



# Proprietary Sample Pretreatment to Detect Living Organisms



What is the relevance of measuring the different ATP portions?


- Extracellular ATP or free ATP is the portion of ATP released by dead cells
- Total ATP is obtained after lysis of the biomass by sonication of the sample
- Intracellular ATP is the portion of ATP from the metabolically active (living) organisms

**Intracellular ATP = Total ATP – Extracellular ATP**

**Bonus Material:** Free ATP can be used in reverse osmosis (RO) applications to monitor contaminant removal

# ATP chemiluminescence reaction

ATP assays using luciferin/luciferase reactions allow to assess microbial load in water sources. The ASTM D4012 (Standard Test Method for Adenosine Triphosphate Content of Microorganisms in Water) was developed as a quick and sensitive alternative to plate counting.



# EZ7300 Series - General

## What?

The first microbiology analyzer using the ATP firefly assay and complying with international standard method **ASTM D4012-81**

The on-line microbiology analyzer brings operators

- Fast results (few minutes  $\approx$  real-time data)
- Automation
- Analytical performance
- Reliability



# Hach EZ7300 Series ATP Analyzer

The first microbiology analyzer using the ATP firefly assay and complying with international standard method **ASTM D4012-81**.



- Fast results (few minutes ≈ real-time data)
  - 10 minutes can be delayed to 2.5 hours
- Analytical performance, LOD of 0.05 pg/mL
- Up to 8 sample streams possible
- Highly sensitive and stable ATP reagents
- Automation and cleaning
  - Rinsing after each sample.
  - Cleaning with HCl to break down any biofilm.
  - Cleaning with NaOH to hydrolyze remaining bacteria

## EZ7300 Series - High analytical performance

- Complete ATP recovery: detection of intracellular, extracellular and total ATP
- Rapid measurement: 7 - 10 minutes (including lysis)
- No delay between sample take-off and measurement
- Smart features: automatic calibration and 3-step cleaning protocol
- Low limit of detection (LOD): 0.05 pg/mL (0.1 pM) ATP \*
- Extended reagent stability (see further)
- Factory configured, tested and calibrated

\* 0.05 pg ~ 50 *E. coli* sized bacteria

## EZ7300 Series – Proprietary sample pretreatment

**Question:** what is one of the main challenges in measuring the microbial load (contamination) in water?

**Answer:** in order to have a clear picture of the microbial load it is important to differentiate ATP portions within living cells from non-living cells


ATP in the water source can be located...

Either inside bacteria or other cells

= intracellular

Or freely in the water surrounding the cells

= extracellular or free ATP





## EZ7300 Series – Proprietary sample pretreatment

### What is the relevance of measuring the different ATP portions?

- Extracellular ATP or free ATP is the portion of ATP released by dead cells
- Total ATP is obtained after lysis of the biomass by sonication of the sample
- Intracellular AP is the portion of ATP from the metabolically active (living) organisms

Or:

$$\text{Intracellular ATP} = \text{Total ATP} - \text{Extracellular ATP}$$

## EZ7300 Series – Proprietary sample pretreatment

**Solution: EZ Series** proprietary sample pretreatment unit performs a controlled ultrasonic lysis of the sample which disrupts the living cells

- Sampling
- Addition luciferase and luciferin
- 1st light output correlates with **free ATP**
- Cell lysis with ultrasonic treatment
- Addition luciferase and luciferin
- 2nd light output correlates with **total ATP**



# Real time monitoring of raw water intake and treated water

The objectives:

## At the Alfred Merritt Smith Water Treatment Facility

- Obtain real-time information to further optimize the ozonation of the raw water coming from Lake Mead.
- Refine and optimize the Chlorine gas disinfection process used on the finished drinking water.



# Saint Paul Regional Water Services' McCarron's Water Treatment Plant

## Water Use (In Million Gallons)

Average daily - (2016) 38.2

## Area Supplied (In Square Miles)

Saint Paul - 56.2

Suburbs - 66.6

## Water Treatment Plant (In Million Gallons)

Daily plant capacity – 144

Total treated, finished water - 88


Number of people served - 432,911




# Optimization of a biological filter and plant risk mitigation

## The objectives:

- Monitor in real time the microbial load of the influent and effluent of their biological filter.
- Determine the efficiency of the biological filter in removing biomass (microbial load)
- Trouble shoot/mitigate risk (corrosion) in other parts of the plant using the grab sample line of the analyzer.



# Field data: drinking water



*The graph shown on the right is the statistical summary from analysis between June and September 2017 regarding filter influent and effluent of a biological filter.*

*Filtration process is effective in removing ~57% of influent biomass based on the ATP data.*


Data courtesy of Saint Paul Regional Water Services, MN, USA

# Drinking water production company senj, vodovod južni ogranak senj

- Water source for production of drinking water is accumulation from two rivers (Lika and Gacka).
- HE Senj - electricity power plant - also use this accumulation for electricity production. Maintenance plan - 3 times a year.
- Maintenance plan affect on quality of raw water, especially on turbidity and microbiology



# Solution



- In normal condition customer use sand filter
- In case of microbiology excess customer will use ATP analyzer as early alarm and switch purification to ultrafiltration



**THANK YOU FOR YOUR TIME TODAY**

